

DISCOVER HOW TO BENEFIT FROM EVERNOA FENO ANALYSER!



Inhale



Exhale in Adult (10s)
or Child (6s) mode



Get the result in 5
seconds

**WE ADAPT TO THE
NEEDS OF YOUR
CLINIC**



Thanks to the system based on customizable KITS*; Select the number of measurements and consume them in **18 months** after the installation of the cassette

*Contact us for more information about the available KITS

**Optimise your time,
consecutive attempts in
<10 seconds!**



Training mode is also
available

The training mode and the measurement of ambient NO do not consume additional measurements

TECHNICAL SPECIFICATIONS



| | | | |
|------------------------|-------------------------------------|--------------------------------------|---|
| Measuring range | 5-300 ppb | Precision | < 2 ppb for values < 50 ppb ≤ 3% for values ≥ 50 ppb |
| Response time | 5 sec | Trueness | ± 4 ppb for values < 50 ppb ≤ 8% for values ≥ 50 ppb |
| Warm-up time | Max. 1 min | Inhalation through the device | Not necessary |
| Portability | Yes | Training | Available |
| Autonomy | > 4h operating > 72h in stand-by | Data storage | 15,000 measures |
| Lifetime | 15,000 measures or 6 years | Ambient NO | Available |
| Size | 240x200x165mm | Weight | 1,4 kg |

And now manage the evolution of your patients with **evernoa view**



Consult with us the possibilities of integration with Electronic Health Record and the management of your database.

IMPORTANT INFORMATION REGARDING EVERNOA

Evernoa measures the fractional concentration of nitric oxide (FENO) from human breath according to the recommendations of the American Thoracic Society (ATS) and the European Respiratory Society (ERS). FENO is a quantitative measure of nitric oxide from the airways that can be used in the asthma diagnosis as well as in its follow-up. Evernoa must be operated by health professionals who have received training and have read the instructions for use. Evernoa can determine FENO in adults and children over 4 years. The operator must differentiate between the measurement mode for children (6s) and the adult mode (10s). Results offered by evernoa should be interpreted by a doctor.



EVERSENS SL
C/ Tajonar 22
Centro Jerónimo de Ayanz, L4
31006, Pamplona-Iruña,
Navarra (Spain)

info@evernoa.com
(+34) 948 16 62 50

